

## 17.1 Caging Groups and Their Photolysis

### Overview of Caged Probes

Flash photolysis of photoactivatable or “caged” probes provides a means of controlling the release — both spatially and temporally — of biologically active products or other reagents of interest.<sup>1–6</sup> The chemical caging process may also confer membrane permeability on the caged ligand, as is the case for caged cAMP<sup>7</sup> and caged luciferin.<sup>8</sup> Molecular Probes’ extensive selection of caged nucleotides, chelators, second messengers and neurotransmitters has tremendous potential for use with live cells.

In addition to caged versions of biologically active molecules, we prepare caged fluorescent dyes that are essentially nonfluorescent until after photolysis. These caged fluorophores have proven useful for photoactivation of fluorescence (PAF) experiments, which are analogous to fluorescence recovery after photobleaching (FRAP) experiments except that the fluorophore is activated upon illumination rather than bleached. Measuring the bright fluorescent signal of the photoactivated fluorophore against a dark background is intrinsically more sensitive than measuring a dark photobleached region against a bright field.<sup>1,9</sup> We have also introduced antibody and streptavidin conjugates of CMNB-caged fluorescein (Section 7.3, Section 7.6) that develop strong green fluorescence following brief illumination with UV light (Figure 7.65).

### Properties of Different Caging Groups

The caging moiety (Table 17.1) is designed to *maximally* interfere with the binding or activity of the molecule. It is detached in microseconds to milliseconds by flash photolysis at  $\leq 360$  nm, resulting in a pulse of the active product. Uncaging can easily be accomplished with UV illumination in the fluorescence microscope or with a UV laser or UV flashlamp. The availability of low-cost UV lasers<sup>10</sup> and UV flashlamps<sup>11</sup> should facilitate photolysis experiments in many laboratories (Table 17.2). A high-resolution, confocal laser-scanning microscope and flash photolysis system for physiological studies has been described.<sup>12</sup> The effects of photolytic release are frequently monitored either with fluorescent probes that measure calcium, pH, other ions or membrane potential, or with electrophysiological techniques.

To date, most of the caged reagents described in the literature have been derivatives of *o*-nitrobenzyl compounds. The nitrobenzyl group is synthetically incorporated into the biologically active molecule by linkage to a hetero-atom (usually *O*, *S* or *N*) as an ether, thioether, ester (including phosphate or thiophosphate esters), amine or similar functional group. Both the structure of the nitrobenzyl compound and the atom to which it is attached affect the efficiency and wavelength required for uncaging. We currently use five different photolabile protecting groups in our caged probes.<sup>2</sup> Their properties are summarized in Table 17.1.

**Table 17.1** Properties of five different caging groups.






Probe	Uncaging Rate *	Photolysis Quantum Yield *	Inertness of Photolysis By-product	Confers Water Solubility	Long-Wavelength Absorption ( $\geq 360$ nm)
CNB	++++	+++++	+++++	+++++	++
NPE	+++	+++	+++	+	++
DMNB	+++	+++	++	+	+++++
DMNPE	+++	+++	+++	+	+++++
CMNB	+++	+++	+	++++	+++

+++++ = Optimal response. + = Poor response.

\* Both the structure of the nitrobenzyl moiety and the atom to which it is attached have some effect on the efficiency and wavelength required for uncaging.

### TECHNICAL NOTE

#### Structures of Five Different Caging Groups



### TECHNICAL NOTE

#### Why Use Caged Probes?

Caged probes permit the rapid release of *natural* biologically active molecules or effectors, (e.g., nucleotides,  $\text{Ca}^{2+}$ , ionophores, inositol 1,4,5-triphosphate and neurotransmitters). When combined with our other probes for signal transduction, caging technology provides results that cannot be obtained by any other means, including patch clamping. It is also possible to “turn on” fluorescence and to create a hapten using caged fluorescent dyes.



## Technical Assistance at Our Web Site ([www.probes.com](http://www.probes.com))

Molecular Probes' Web site contains the complete text of this Handbook as well as a number of other features to assist the researcher, including:

- **Product searches** by product name or catalog number
- **Color photomicrographs** that show our products in action
- **Bibliographies** for all products for which we have references
- **Keyword searches** of our entire bibliography of over 44,000 references
- **Product information sheets** for many kits and reagents
- **Technical bulletins**, including *BioProbes* newsletters and other product literature
- **Chemical structures, technical data and material safety data sheets**

Check our Web site frequently for information on our newest products and the most recent additions to our bibliography, as well special offers on featured products.

Additional information on the scientific and technical background of our products is available from our Technical Assistance Department:

### In the U.S. and Canada

Phone: 541-465-8353

Fax: 541-465-4593

E-mail: [tech@probes.com](mailto:tech@probes.com)

### In Europe

Phone: +31-71-5233431

Fax: +31-71-5241883

E-mail: [eurotech@probes.nl](mailto:eurotech@probes.nl)

- Probes caged with our patented  $\alpha$ -carboxy-2-nitrobenzyl (CNB) caging group<sup>13</sup> generally have the most advantageous properties. These include good water solubility, very fast uncaging rates in the microsecond range, high photolysis quantum yields (from 0.2–0.4) and biologically inert photolytic by-products. Although the absorption maximum of the CNB caging group is near 260 nm, its absorption spectrum tails out to approximately 360 nm, allowing successful photolysis using light with wavelengths  $\leq 360$  nm.
- The 1-(2-nitrophenyl)ethyl (NPE) caging groups have properties similar to those of CNB and can also be photolyzed at  $\leq 360$  nm.
- As compared with CNB and NPE, the 4,5-dimethoxy-2-nitrobenzyl (DMNB) and 1-(4,5-dimethoxy-2-nitrophenyl)ethyl (DMNPE) caging groups have longer-wavelength absorption (absorption maximum  $\sim 355$  nm) and therefore absorb 340–360 nm light more efficiently. However, photolysis rates and quantum yields of DMNB- and DMNPE-caged probes are generally lower than those obtained for CNB-caged probes.
- The 5-carboxymethoxy-2-nitrobenzyl (CMNB) caging group provides an absorption maximum of intermediate wavelength (absorption maximum  $\sim 310$  nm), while imparting significant water solubility to the caged probe. Its photolysis rate and quantum yield are intermediate between those of CNB- and DMNB-caged probes.

Experiments utilizing probes caged with any of the above caging groups, except the CNB caging group, may require the addition of dithiothreitol (DTT, D-1532; Section 2.1). This reducing reagent prevents the potentially cytotoxic reaction between amines and the 2-nitrosobenzoyl photolytic by-products.<sup>14</sup>

## References

1. *Methods Enzymol* 291, 63 (1998);
2. *Methods Enzymol* 291, 30 (1998);
3. *Curr Opin Neurobiol* 6, 379 (1996);
4. *Biological Applications of Photochemical Switches*, Morrison H, Ed. pp. 243–305 (1993);
5. *Optical Microscopy: Emerging Methods and Applications*, Herman B, Lemasters JJ, Eds. pp. 27–85 (1993);
6. *Annu Rev Physiol* 55, 755 (1993);
7. *Nature* 310, 74 (1984);
8. *Biotechniques* 15, 848 (1993);
9. *Cell Biology: A Laboratory Handbook*, 2nd Ed., Vol. 3, Celis JE, Ed. pp. 127–135 (1998);
10. *J Neurosci Methods* 66, 47 (1996);
11. *Pflugers Arch* 411, 200 (1988);
12. *Cell Calcium* 21, 441 (1997);
13. *US* 5,635,608;
14. *Annu Rev Biophys Biophys Chem* 18, 239 (1989).

**Table 17.2** Suppliers of instrumentation for photolysis of caged compounds.

Company	Location	Web Site
Cairn Research Ltd.	Faversham, UK	<a href="http://www.cairnweb.com">www.cairnweb.com</a>
Intracellular Imaging, Inc.	Cincinnati, OH, USA	<a href="http://www.intracellular.com">www.intracellular.com</a>
Fryer Company, Inc.	Huntley, IL, USA	<a href="http://www.fryerco.com/prairie">www.fryerco.com/prairie</a>
Hi-Tech Scientific	Salisbury, UK	<a href="http://www.hi-techsci.co.uk">www.hi-techsci.co.uk</a>
Photonic Instruments, Inc.	Arlington Heights, IL, USA	<a href="http://www.photonic-instruments.com">www.photonic-instruments.com</a>
Prairie Technologies, LLC	Middleton, WI, USA	<a href="http://www.prairie-technologies.com">www.prairie-technologies.com</a>
Rapp OptoElectronic	Hamburg, Germany	<a href="http://www.rapp-opto.com">www.rapp-opto.com</a>
T.I.L.L. Photonics	Martinsried, Germany	<a href="http://www.till-photonics.com">www.till-photonics.com</a>

## Extra Copies of this Handbook

Publication of this *Handbook* involves considerable effort and expense on the part of Molecular Probes. Single copies are provided free upon request to new customers. Additional copies can be ordered from Molecular Probes under catalog number H-19999. In Japan, Korea, Australia and Mexico, copies of the *Handbook* must be requested from one of our distributors in those countries (see page iv). Copies may also be requested from our other distributors (see page iv).