

cal arm of the immune system, as well as the evolution of antibodies and the role they may play in human diseases.^{9,10} Several reviews have discussed the chemistry of the different reactive oxygen species and their lipid peroxidation products.^{11–16}

The importance of the nitric oxide radical (abbreviated NO) and other reactive oxygen species as biological messengers has been increasingly recognized during the last several years.^{17–22} Section 19.3 is devoted to our probes for promoting, inhibiting or detecting nitric oxide production in a variety of experimental systems.

References

1. *Prog Neurobiol* 62, 633 (2000);
2. *Annu Rev Biochem* 64, 97 (1995);
3. *Methods Enzymol* 234, (Complete Volume) (1994);
4. *Methods Enzymol* 233, (Complete Volume) (1994);
5. *Methods Enzymol* 186, (Complete Volume) (1990);
6. *Methods Enzymol* 105, (Complete Volume) (1984);
7. *J Biol Chem* 268, 18502 (1993);
8. *Pharmacol Ther* 53, 375 (1992);
9. *Science* 293, 1806 (2001);
10. *Proc Natl Acad Sci U S A* 97, 10930 (2000);
11. *Free Radic Biol Med* 30, 463 (2001);
12. *Clin Chem* 41, 1819 (1995);
13. *Free Radic Biol Med* 19, 77 (1995);
14. *Free Radic Biol Med* 18, 1033 (1995);
15. *Free Radic Biol Med* 18, 775 (1995);
16. *Biochimie* 76, 355 (1994);
17. *Curr Biol* 7, R376 (1997);
18. *J Med Chem* 38, 4343 (1995);
19. *Annu Rev Biochem* 63, 175 (1994);
20. *Cell* 78, 919 (1994);
21. *Am J Physiol* 262, G379 (1992);
22. *Science* 257, 494 (1992).

19.2 Generating and Detecting Reactive Oxygen Species

Molecular Probes prepares an assortment of probes for the generation of reactive oxygen species (ROS) — singlet oxygen, hydroxyl radicals, superoxide, hydroperoxides and peroxides (Table 19.1) — as well as for their fluorometric detection in solution. Although there are no equilibrium sensors that continuously monitor the level of reactive oxygen species, this section discusses a number of probes that trap or otherwise react with singlet oxygen, hydroxyl radicals or superoxide. The optical or electron spin properties of the resulting products are in some way a measure of the presence or quantity of the reactive oxygen species and, in some cases, can measure the kinetics and location of their formation.

Generating Singlet Oxygen

Singlet oxygen is responsible for much of the physiological damage caused by reactive oxygen species, including nucleic acid modification through selective reaction with deoxyguanosine.¹ The lifetime of singlet oxygen is sufficiently long (4.4 microseconds in water²) to permit significant diffusion in cells and tissues.³ In the laboratory, singlet oxygen is usually generated in one of three ways: photochemically from dioxygen (³O₂) using a

photosensitizing dye;⁴ chemically, either by thermal decomposition of a peroxide or diacetate; or by microwave discharge through an oxygen stream. Singlet oxygen can be detected by its characteristic weak chemiluminescence at 1268 nm⁵ or at 634 and 703 nm.⁶

Hypocrellins and Hypericin

Among the most efficient reagents for generating singlet oxygen are the photosensitizers hypocrellin A (H-7515, Figure 18.10), hypocrellin B (H-7516, Figure 18.11) and hypericin (H-7476, Figure 18.9). These heat-stable dyes exhibit quantum yields for singlet oxygen generation in excess of 0.7, as well as high photostability, making them important agents for both anti-cancer and antiviral therapy.^{7–13} Hypocrellins are also specific and potent inhibitors of protein kinase C^{14,15} (Section 18.3), whereas hypericin is an effective inhibitor of both protein kinase C and tyrosine protein kinase¹⁶ with antiretroviral activity.¹⁷

Because their chemical reactivities are well characterized,^{18,19} hypocrellins and hypericin are amenable to conjugation to a variety of primary and secondary detection reagents. Not only do these photosensitizing dyes efficiently oxidize diaminobenzidine (DAB) to form an insoluble, electron-dense DAB oxidation product (see




Figure 19.1 R-14000 rose bengal diacetate.




Figure 19.2 Reaction of *trans*-1-(2'-methoxyvinyl)pyrene (M-7913) with singlet oxygen (¹O₂), yielding a diacetate intermediate that generates chemiluminescence (CL) upon decomposition to 1-pyrenecarboxaldehyde.

Fluorescent Probes for Photoconversion of Diaminobenzidine Reagents in Section 1.5), but they exhibit modest fluorescence quantum yields and broad UV and visible spectra. Thus, hypocrelin- and hypericin-labeled detection reagents are compatible with fluorescence, light or electron microscopy applications.²⁰

Halogenated Fluoresceins and Rhodamines

Halogenated derivatives of fluorescein dyes are known to be effective photosensitizers and singlet oxygen generators.²¹ Bioconjugates of eosin and erythrosin fluorophores can be prepared using the reactive derivatives described in Chapter 1, Chapter 2 and Chapter 3. The practical use of eosin and other conjugates to improve ultrastructural resolution through singlet oxygen generation is described in Fluorescent Probes for Photoconversion of Diaminobenzidine Reagents (see Section 1.5). The generation of singlet oxygen can be targeted to the mitochondria by several cationic dyes,²² such as tetrabromorhodamine 123 (T-7539), which has a quantum yield for singlet oxygen generation of 0.65–0.7.^{22,23}

Rose Bengal Diacetate

Rose bengal diacetate (R-14000) is an efficient, cell-permeant generator of singlet oxygen.²⁴ It is an iodinated xanthene derivative that has been chemically modified by the introduction of acetate groups (Figure 19.1). These modifications inactivate both its fluorescence and photosensitization properties, while increasing its ability to cross cell membranes. Once inside a live cell, esterases remove the acetate groups and restore rose bengal to its native structure. This intracellular localization allows rose bengal diacetate to be a very effective photosensitizer.

Merocyanine 540

Photolysis of merocyanine 540 (M-24571) produces both singlet oxygen and other reactive oxygen species, including oxygen radicals.^{25–28} Merocyanine 540 is commonly used as a photosensitizer in photodynamic therapy.^{10,29–37}

Detecting Singlet Oxygen

trans-1-(2'-Methoxyvinyl)pyrene (M-7913) can be used to detect picomole quantities of singlet oxygen in chemical and biological systems (Figure 19.2), making this compound perhaps the most sensitive singlet oxygen probe currently available.^{38–40} Furthermore, this highly selective chemiluminescent probe does not react with other activated oxygen species such as hydroxyl radical, superoxide or hydrogen peroxide.

Generating Hydroxyl and Superoxide Radicals

Hydroxyl and superoxide radicals have been implicated in a number of pathological conditions, including ischemia, reperfusion and aging. The superoxide anion (Table 19.1) may also play a role in regulating normal vascular function. The hydroxyl radical is a very reactive oxygen species that has a lifetime of about 2 nanoseconds in aqueous solution and a radius of diffusion of about 20 Å. Thus, it induces peroxidation only when it is generated in close proximity to its target. The hydroxyl radical can be derived from superoxide in a reaction catalyzed by Fe²⁺ or other transition metals, as well as by the effect of ionizing radiation on dioxygen. Superoxide is most effectively generated from a hypo-

xanthine–xanthine oxidase generating system.^{41–43} In one cell-based assay for hydroxyl radical formation, D-phenylalanine is specifically converted to D-tyrosine.⁴⁴ In phagocytic cells, H₂O₂ also produces *O,O'*-dityrosine, an oxidative crosslink product of appropriately situated tyrosine residues, which is formed through the intermediacy of phenoxyl radicals.⁴⁵

Malachite Green

Malachite green is a nonfluorescent photosensitizer that absorbs at long wavelengths (~630 nm). Its photosensitizing action can be targeted to particular cellular sites by conjugating malachite green isothiocyanate (M-689, Figure 1.83) to specific antibodies.^{46,47} Enzymes and other proteins within ~10 Å of the binding site of the malachite green-labeled antibody can then be selectively destroyed upon irradiation with long-wavelength light. Studies by Jay and colleagues have demonstrated that this photo-induced destruction of enzymes in the immediate vicinity of the chromophore is apparently the result of localized production of hydroxyl radicals, which have short lifetimes that limit their diffusion from the site of their generation.⁴⁸

1,10-Phenanthroline Iodoacetamide

Conjugation of the iodoacetamide of 1,10-phenanthroline (P-6879) to thiol-containing ligands confers the metal-binding properties of this important complexing agent on the ligand. For example, the covalent copper–phenanthroline complex of oligonucleotides or nucleic acid-binding molecules in combination with hydrogen peroxide acts as a chemical nuclease to selectively cleave DNA or RNA.^{49–54} Hydroxyl radicals or other reactive oxygen species appear to be involved in this cleavage.^{55,56}

Detecting Hydroxyl and Superoxide Radicals

Fluorogenic Spin Traps

Hydroxyl radicals have usually been detected after reaction with spin traps. We offer TEMPO-9-AC (A-7923) and proxyl fluorescamine^{57–60} (C-7924), two fluorogenic probes for detecting hydroxyl radicals and superoxide. Each of these molecules contains a nitroxide moiety that effectively quenches its fluorescence. However, once TEMPO-9-AC or proxyl fluorescamine traps a hydroxyl radical or superoxide, its fluorescence is restored and the radical's electron spin resonance signal is destroyed, making these probes useful for detecting radicals either by fluorescence or by electron spin resonance spectroscopy. Proxyl fluorescamine can also be used to detect the methyl radicals that are formed by the reaction of hydroxyl radicals with DMSO.⁶⁰ Radical-specific scavengers — such as the superoxide-specific *p*-benzoquinone and superoxide dismutase⁶¹ or the hydroxyl radical-specific mannitol and dimethylsulfoxide (DMSO)^{57,62,63} — can be used to identify the detected species.

Chemiluminescent and Chromogenic Reagents for Detecting Superoxide

In the absence of apoaquorin, the luminophore coelenterazine (C-2944) produces chemiluminescence in response to superoxide generation in phorbol ester- or chemotactic peptide-stimulated neutrophils.⁶⁴ Unlike luminol, coelenterazine exhibits luminescence that does not depend on the activity of cell-derived myeloperoxidase and is not inhibited by azide.⁶⁴

In addition to coelenterazine, we offer MCLA (M-23800) for detecting superoxide. MCLA and coelenterazine are superior alternatives to lucigenin (L-6868) for this application because lucigenin can reportedly sensitize superoxide production, leading to false-positive results.^{65–69} An additional advantage of MCLA is that its pH optimum for luminescence generation is closer to the physiological near-neutral range than are the pH optima of luminol and lucigenin.⁷⁰

Lucigenin (L-6868) exhibits chemiluminescence that is reported to be sensitive to the superoxide anion.^{71–75} Lucigenin has been employed to investigate superoxide generation in spermatozoa,⁷⁶ L929 cells,⁷⁷ chondrocytes contained within the matrix of living cartilage tissue,⁷⁸ and mitochondria of alveolar macrophages.⁷⁵ It has also been used to examine the role of the superoxide anion in reoxygenation injury in isolated rat hepatocytes.^{79,80} We have purified our lucigenin to remove a blue-fluorescent impurity that is found in some commercial samples.

Nitro blue tetrazolium salt (NBT, N-6495; Table 19.2) and other tetrazolium salts are chromogenic probes useful for superoxide determination.^{81,82} These probes are also widely used for detecting redox potential of cells for viability, proliferation and cytotoxicity assays; see below for more information.

Detecting Hydrogen Peroxide, Hydroperoxides and Peroxyl Radicals

In peroxisomes, H₂O₂ is produced by several enzymes that use molecular oxygen to oxidize organic compounds. This H₂O₂ is then utilized by catalase to oxidize other substrates, including

phenols, formic acid, formaldehyde and alcohol. In liver and kidney cells, these oxidation reactions are important for detoxifying a variety of compounds in the bloodstream.^{83–86} However, H₂O₂ also plays a role in neurodegenerative and other disorders through induction of apoptosis,⁸⁷ oxidation of glutathione, modification of intracellular Ca²⁺ levels and mitochondrial potential and induction of DNA strand breaks.⁸⁸ In addition, H₂O₂ is released from cells during hypoxia.⁸⁹

Peroxidation of unsaturated lipids plays an important role in cell membrane properties,⁹⁰ signal transduction pathways,^{91,92} apoptosis and the deterioration of foods and other biological compounds.⁹³ Lipid oxidation may also be responsible for aging, as well as pathological processes such as drug-induced phototoxicity and atherosclerosis.⁹⁴ Lipid hydroperoxides have been reported to accumulate in oxidatively stressed individuals, including HIV-infected patients.⁹⁵ To directly assess the extent of lipid peroxidation, researchers either measure the amount of lipid hydroperoxides directly or detect the presence of secondary reaction products^{96–98} (e.g., 4-hydroxy-2-nonenal or malonaldehyde; see below).

Peroxyl radicals are formed by the decomposition of hydroperoxides, including lipid hydroperoxides. Experimentally, peroxyl radicals are generated from hydroperoxides such as cumene hydroperoxide or from compounds such as 2,2'-azobis(2-amidinopropane). Free radical-mediated damage in cells is often the result of lipid peroxidation. A problem with investigating the link between lipid peroxidation and diseases such as atherosclerosis, diabetes and Parkinson's has been the lack of suitable methods to detect the relationship between lipid peroxidation and the onset of such diseases.⁹⁹

Table 19.2 Tetrazolium salts for detecting redox potential in living cells and tissues.

| Cat # | Tetrazolium Salt | Color of Formazan | Water Solubility of Formazan | Applications |
|--------------|---|-------------------|------------------------------|--|
| M-6494 (MTT) | 3-(4,5-Dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide | purple | no | <ul style="list-style-type: none"> • Superoxide generation by fumarate reductase¹ and nitric oxide synthase² • Mitochondrial dehydrogenase activity³ • Cell viability and proliferation^{4–9} • Neuronal cell death¹⁰ • Platelet activation¹¹ • Tumor cell adhesion¹² and invasion¹³ • Multidrug resistance¹⁴ • <i>In vitro</i> toxicity testing^{15–17} |
| N-6495 (NBT) | Nitro blue tetrazolium chloride | deep blue | no | <ul style="list-style-type: none"> • Superoxide generation by xanthine oxidase¹⁸ • Neutrophil oxidative metabolism^{19,20} • NADPH diaphorase activity^{21–23} • Succinic dehydrogenase histochemistry²⁴ |
| X-6493 (XTT) | 2,3-Bis-(2-methoxy-4-nitro-5-sulfophenyl)-2 <i>H</i> -tetrazolium-5-carboxanilide | orange | yes | <ul style="list-style-type: none"> • Antifungal susceptibility²⁵ • Drug sensitivity of cells²⁶ • Parasitic nematode viability²⁷ • Tumor cell cytotoxicity²⁸ |

1. J Biol Chem 270, 19767 (1995); 2. J Biol Chem 269, 12589 (1994); 3. Cytometry 13, 532 (1992); 4. Biotechniques 25, 622 (1998); 5. J Immunol Methods 168, 253 (1994); 6. Anal Biochem 214, 190 (1993); 7. J Immunol Methods 164, 149 (1993); 8. Anal Biochem 205, 8 (1992); 9. J Immunol Methods 89, 271 (1986); 10. J Cell Biol 128, 201 (1995); 11. J Immunol Methods 159, 253 (1993); 12. J Immunol Methods 164, 255 (1993); 13. Cancer Res 54, 3620 (1994); 14. Leuk Res 16, 1165 (1992); 15. Biosci Biotechnol Biochem 56, 1472 (1992); 16. J Immunol Methods 144, 141 (1991); 17. J Immunol Methods 131, 165 (1990); 18. J Reprod Fertil 97, 441 (1993); 19. Clin Chim Acta 221, 197 (1993); 20. J Leukoc Biol 53, 404 (1993); 21. Neurosci Lett 155, 61 (1993); 22. Proc Natl Acad Sci U S A 88, 7797 (1991); 23. Proc Natl Acad Sci U S A 88, 2811 (1991); 24. Histochemistry 76, 381 (1982); 25. Antimicrob Agents Chemother 36, 1619 (1992); 26. Cancer Res 48, 4827 (1988); 27. Parasitology 107, 175 (1993); 28. J Immunol Methods 147, 153 (1992).

Amplex Red Reagent

In the presence of horseradish peroxidase (HRP), the Amplex Red reagent (10-acetyl-3,7-dihydroxyphenoxazine, A-12222, A-22177) reacts with H_2O_2 with a 1:1 stoichiometry to produce highly fluorescent resorufin (Figure 10.50). The Amplex Red reagent has greater stability, yields less background and produces a red-fluorescent product that is more readily detected than the similar reduced methylene blue derivatives commonly used for colorimetric determination of lipid peroxides in plasma, sera, cell extracts and a variety of membrane systems.^{33,100,101} Using the Amplex Red reagent in conjunction with HRP, we have found that release of hydrogen peroxide to the medium by as few as 2000 phorbol ester-stimulated neutrophils can be detected in a fluorescence microplate reader. The Amplex Red reagent has been used to detect H_2O_2 release from human keratinocytes after UV-B radiation¹⁰² and from human leukocytes after phorbol myristate acetate treatment¹⁰³ and from proliferating U937 cells.¹⁰⁴ Using the Amplex Red reagent, researchers have discovered that antibodies can convert molecular oxygen to H_2O_2 , which may be important in understanding a new chemical arm of the immune system, as well as the evolution of antibodies and the role they may play in human diseases.^{105,106}

Amplex Red Hydrogen Peroxide/Peroxidase Assay Kit

Our Amplex Red Hydrogen Peroxide/Peroxidase Assay Kit (A-22188), which is described in greater detail in Section 10.5, includes the reagents required for analysis of either H_2O_2 in solutions or peroxidase activity in an ultrasensitive, one-step assay that can be measured using either fluorescence or absorbance. The Amplex Red peroxidase substrate can detect the presence of active peroxidases and the release of H_2O_2 from biological samples, including cells and cell extracts^{103,107–109} and is also useful for detecting H_2O_2 that is produced as a product of enzyme-coupled reactions.^{110,111} The Amplex Red Hydrogen Peroxide/Peroxidase Assay Kit contains:

- The Amplex Red reagent
- Horseradish peroxidase (HRP)
- H_2O_2 for use as a positive control
- Concentrated reaction buffer
- A detailed protocol

Each kit provides a sufficient amount of reagent for approximately 500 assays using a reaction volume of 100 μL per assay.

Several additional kits that utilize our Amplex Red peroxidase substrate to detect H_2O_2 in coupled enzymatic reactions are described in Section 10.5.

Amplex Red Xanthine/Xanthine Oxidase Assay Kit

Xanthine oxidase plays a key role in the production of free radicals, including superoxide in the body. In healthy individuals, xanthine oxidase is present in appreciable amounts only in the liver and the jejunum. However, in various liver disorders the enzyme is released into circulation. Therefore, determination of serum xanthine oxidase level serves as a sensitive indicator of acute liver damage such as jaundice.¹¹² Previously, researchers have used chemiluminescence^{42,113} or absorbance¹¹⁴ to monitor xanthine oxidase activity. Using the Amplex Red Xanthine/Xanthine Oxidase Assay Kit (A-22182), one can detect xanthine oxidase activity in a purified system at levels as low as 0.1 mU/mL by fluorescence (Figure 19.3). The kit can also detect as little as 200 nM hypoxanthine or xanthine (Figure 19.4).

Xanthine oxidase catalyzes the oxidation of hypoxanthine or xanthine to uric acid and superoxide. In the reaction mixture, the superoxide spontaneously degrades to hydrogen peroxide, which in the presence of horseradish peroxidase reacts stoichiometrically with the Amplex Red reagent to generate the red-fluorescent oxidation product resorufin. Resorufin has absorption and emission maxima at approximately 563 and 587 nm, respectively (Figure 10.5), and because the extinction coefficient is high ($54,000 \text{ cm}^{-1}\text{M}^{-1}$), the assay can be performed either fluorometrically or spectrophotometrically.

The Amplex Red Xanthine/Xanthine Oxidase Assay Kit (A-22182) contains:

- The Amplex Red reagent
- Dimethylsulfoxide (solvent for the Amplex Red reagent)
- Horseradish peroxidase (HRP)
- Hydrogen peroxide
- Concentrated reaction buffer
- Xanthine oxidase from buttermilk
- A hypoxanthine solution
- A xanthine solution
- A detailed protocol




Figure 19.3 Detection of xanthine oxidase using the Amplex Red reagent-based assay (A-22182). Each reaction contained 50 μM Amplex Red reagent, 0.2 U/mL horseradish peroxidase, 0.1 mM hypoxanthine and the indicated amount of xanthine oxidase in 1 \times reaction buffer. After 30 minutes, fluorescence was measured in a fluorescence microplate reader using excitation at $530 \pm 12.5 \text{ nm}$ and detection at $590 \pm 17.5 \text{ nm}$. A background of 65 fluorescence units was subtracted from each data point. The inset shows the assay's sensitivity and linearity at low enzyme concentrations.




Figure 19.4 Detection of hypoxanthine using the Amplex Red reagent-based assay (A-22182). Each reaction contained 50 μM Amplex Red reagent, 0.2 U/mL horseradish peroxidase, 20 mU/mL xanthine oxidase and the indicated amount of hypoxanthine in 1 \times reaction buffer. Reactions were incubated at 37°C. After 30 minutes, fluorescence was measured in a fluorescence microplate reader using excitation at $530 \pm 12.5 \text{ nm}$ and detection at $590 \pm 17.5 \text{ nm}$. A background of 54 fluorescence units was subtracted from each data point. The inset shows the assay's sensitivity and linearity at low enzyme concentrations.

Amplex Red reagent-based assays for numerous enzymes and analytes are described in Section 10.5. A particularly important advantage of these assays is their use of natural substrates, rather than synthetic pseudosubstrates, in the assays.

Each kit provides sufficient reagents for ~400 assays using either a fluorescence or absorbance microplate reader and reaction volumes of 100 μ L per assay.

***cis*-Parinaric Acid**

Fluorescence quenching of the fatty acid analog *cis*-parinaric acid (P-1901) has been used in several lipid peroxidation assays,^{99,115–119} including quantitative determinations in living cells.^{120–125} In a study investigating the membrane antioxidant properties of the *bcl-2* proto-oncogene product, researchers used *cis*-parinaric acid to detect lipid hydroperoxides together with 6-carboxy-2',7'-dichlorodihydrofluorescein diacetate, di(acetoxymethyl ester) (C-2938) to detect cytosolic reactive oxygen molecules.¹²⁶ Parinaric acid's extensive unsaturation (Figure 13.11) makes it quite susceptible to oxidation if not rigorously protected from air.¹²⁷ Consequently, we offer *cis*-parinaric acid specially packaged in 10 ampules sealed under argon, each containing 10 mg (P-1901). During experiments, we advise handling parinaric acid samples under inert gas and preparing solutions using degassed buffers and solvents. Parinaric acid is also somewhat photolabile and undergoes photodimerization when exposed to intense illumination, resulting in loss of fluorescence.¹²⁸

Diphenyl-1-Pyrenylphosphine

The direct fluorometric detection of hydroperoxides in lipids, serum, tissues and foodstuffs^{129,130} is now possible using the reagent diphenyl-1-pyrenylphosphine (DPPP, D-7894). DPPP is essentially nonfluorescent until oxidized to a phosphine oxide by peroxides. DPPP has previously been used to detect picomole levels of hydroperoxides by HPLC.^{131–134} Its solubility in lipids makes DPPP quite useful for detecting hydroperoxides in the membranes of living cells¹³⁵ and in low-density lipoprotein particles,¹³⁶ as well as for studying superoxide dismutase (SOD), which catalyzes the conversion of superoxide to hydrogen peroxide.

BODIPY Dyes: Peroxyl Radical Scavengers

A unique assay for peroxyl radicals uses some of our BODIPY dyes, including the BODIPY 581/591 fatty acid⁹⁰ (D-3861) and the BODIPY 665/676 dye (B-3932) to measure the antioxidant activity in an organic lipid environment or in a liposomal medium.^{90,137,138} This assay is based on the loss or shift of the dye's fluorescence as the result of its interaction with peroxyl radicals, and on the retention of the fluorescent signal in the presence of antioxidants that intercept these free radicals. It has been proposed that this assay is suitable for a fluorescence-based microplate reader or flow cytometer to examine the effect of lipid peroxidation on a cell-by-cell basis.¹³⁸

The oxidation sensitivity of the conjugated double bonds of the BODIPY 581/591 fatty acid (Figure 19.5) has also been exploited in a ratiometric assay of lipid oxidation in live cells that is essentially independent of uneven dye loading, cell thickness and compartmentalization.¹³⁹ In this assay, the cell membranes of rat-1 fibroblasts and myocardial cells were labeled with the BODIPY 581/591 fatty acid, and then cumene hydroperoxide was added to induce oxidation. Upon oxidation, the BODIPY 581/591 fluorophore exhibited a shift in its fluorescence from red to green, which was observed using confocal laser-scanning microscopy. As compared with arachidonic acid, the fluorescence of this BODIPY probe was reportedly twice as sensitive to oxidation.

Peroxy radicals have also been detected in erythrocyte and red blood cell membranes using BODIPY FL EDA¹¹⁷ (D-2390, Section 3.3), a water-soluble BODIPY dye, or BODIPY FL hexadecanoic acid (D-3821, Section 13.2). BODIPY FL hexadecanoic acid exhibits the red shift common to the fluorescence of lipophilic BODIPY dyes when they are concentrated (Figure 13.6), permitting ratiometric measurements of hydroxyl radical production and allowing the onset of lipid peroxidation in live cells to be monitored.¹¹⁷

Other Scavengers for Peroxyl Radicals

The fluorescence of several other probes is lost following interaction with peroxyl radicals. Lipophilic fluorescein dyes such as hexadecanoylamino fluorescein⁹⁹ (H-110, Section 13.5) and fluorescein-labeled phosphatidylethanolamine (F-362, Section 13.2) have been or potentially are useful for detecting peroxyl radical formation in membranes and in solution. Phycobiliproteins such as B-phycoerythrin,^{140–145} R-phycoerythrin^{99,143,146,147} and allophycocyanin¹⁴⁸ (P-800, P-801, A-803, A-819; Section 6.4) may be similarly useful. R-phycoerythrin has been used to detect and measure total plasma antioxidant capacity, including peroxyl radicals.^{143,145,147,149}

Luminol

Although luminol (L-8455) is not useful for detecting superoxide in living cells,⁷³ it is commonly employed to detect peroxidase- or metal ion-mediated oxidative events.^{150–152} Used alone, luminol can detect oxidative events in cells rich in peroxidases, including granulocytes^{153–156} and spermatozoa.⁷⁶ This probe has also been used in conjunction with horseradish peroxidase (HRP) to investigate reoxygenation injury in rat hepatocytes.^{79,157} In these experiments, it is thought that the primary species being detected is hydrogen peroxide. In addition, luminol has been employed to detect peroxynitrite,^{158–160} a molecule thought to be generated in a variety of pathological conditions.⁶² Phospholipid hydroperoxides have been determined directly by chemiluminescence-detected HPLC with luminol^{161–163} (L-8455) or a combination of luminol and cytochrome *c*.¹⁶⁴ This chemiluminescent probe has also been employed as a chemical sensor for dioxygen and nitrogen dioxide. It has been reported that luminol's chemiluminescence response to oxidative species may be competitively inhibited by biomolecules containing sulfhydryl and thioether groups.^{165,166}

Detecting 4-Hydroxy-2-Nonenal

Formation of 4-hydroxy-2-nonenal from linoleic acid is a major cause of lipid peroxidation-induced liver toxicity. Several reagents for the direct fluorometric detection of aldehydes are described in Section 3.2. The modification of 4-hydroxy-2-nonenal or malonaldehyde with a fluorescent or chromophoric hydrazine reagent coupled with the separation and detection of the reaction product by a chromatography-based technique has rarely been reported but appears to be one of the most promising approaches for detecting these lipid peroxidation products.

Assaying Oxidative Activity in Live Cells with Leuco Dyes

Assaying oxidative activity in living cells with fluorogenic, chemiluminescent or chromogenic probes is complicated by the

possibility of having multiple forms of reactive oxygen in the same cell. In addition, the nitric oxide radical (Section 19.3) may produce the same changes in the optical properties of the probe as do other reactive oxygen molecules. Blocking agents and enzyme inhibitors can sometimes help to sort out the species responsible for the probe's optical response. Quantitative analysis is also difficult because of: 1) the high intracellular concentration of glutathione, which can form thiyl or sulfanyl radicals or otherwise trap or reduce oxygen species;⁴¹ 2) the variable concentration of metals, which can either catalyze or inhibit radical reactions; and 3) the presence of other free radical-quenching agents such as spermine.¹⁶⁷

Fluorescein, rhodamine and various other dyes can be chemically reduced to colorless, nonfluorescent leuco dyes. These "dihydro" derivatives are readily oxidized back to the parent dye by some reactive oxygen species and thus can serve as fluorogenic probes for detecting oxidative activity in cells and tissues;^{168–170} however, their oxidation may not easily discriminate between the various reactive oxygen species. It has been reported that dihydroethidium, dichlorodihydrofluorescein (H₂DCF) and dihydro-rhodamine 123 react with intracellular hydrogen peroxide in a reaction mediated by peroxidase, cytochrome *c* or Fe²⁺.^{170–173} The leuco dyes also serve as fluorogenic substrates for peroxidase enzymes (Section 10.5). All of these reagents are slowly oxidized by air back to the parent fluorescent dyes, and in some cases light appears to accelerate their oxidation.

Dichlorodihydrofluorescein Diacetate and Its Analogs

Dichlorodihydrofluorescein diacetate (H₂DCFDA, D-399; Figure 19.6), also known as dichlorofluorescein diacetate, is commonly used to detect the generation of reactive oxygen intermediates in neutrophils and macrophages.^{174–178} Cell-permeant H₂DCFDA may also be extremely useful for assessing the overall oxidative stress in toxicological phenomena.^{169,170} A review by Tsuchiya and colleagues outlined methods for visualizing the generation of oxidative species in whole animals. For example, they suggest using propidium iodide (P-1304, P-3566; FluoroPure Grade, P-21493; Section 8.1) with H₂DCFDA to simultaneously monitor oxidant production and cell injury.¹⁷⁹ H₂DCFDA has been used to visualize oxidative changes in carbon tetrachloride-perfused rat liver¹⁸⁰ and in venular endothelium during neutrophil activation,¹⁸¹ as well as to examine the effect of ischemia and reperfusion in lung and heart tissue.^{182,183} Using H₂DCFDA, researchers characterized hypoxia-dependent peroxide production in *Saccharomyces cerevisiae* as a possible model for ischemic tissue destruction.¹⁸⁴ A variety of toxicological phenomena in cultured cells have also been investigated with H₂DCFDA, including:

- Amyloid β protein-mediated increases in hydrogen peroxide in PC12 cells¹⁸⁵
- Effects of calcium antagonists on oxidative metabolism in dissociated rat cerebellar and cortical neurons¹⁸⁶
- Methamphetamine-induced oxidative stress in dopaminergic neurons¹⁸⁷
- Effect of lipopolysaccharides on the level of oxygen metabolites in rat liver Kupffer cells¹⁸⁸
- Nephrotoxin-induced oxidative stress in isolated proximal tubular cells¹⁸⁹
- Nickel-induced increases in oxidant levels in Chinese hamster ovary (CHO) cells¹⁹⁰
- Effect of transforming growth factor- β 1 on the overall oxidized state of mouse osteoblastic cells¹⁹¹

In neutrophils, H₂DCFDA has proven useful for flow cytometric analysis of nitric oxide, forming a product that has spectral properties identical to those produced when it reacts with hydrogen peroxide.¹⁹² In this study, H₂DCFDA's reaction with nitric oxide was blocked by adding the arginine analog *N*^G-methyl-L-arginine (L-NMMA, M-7898; Section 19.3) to the cell suspension.¹⁹² H₂DCF is not directly oxidized by either superoxide or free hydroxyl radical in aqueous solution.¹⁷⁰ However, it has been reported that 2',7'-dichlorofluorescein — the oxidation product of H₂DCF — can be further oxidized to a phenoxyl radical in a horseradish peroxidase-catalyzed reaction, and this reaction may complicate the interpretation of results obtained with this probe in cells undergoing oxidative stress.¹⁹³




Figure 19.5 D-3861 4,4-difluoro-5-(4-phenyl-1,3-butadienyl)-4-bora-3a,4a-diaza-s-indacene-3-undecanoic acid.




Figure 19.6 D-399 2',7'-dichlorodihydrofluorescein diacetate.

TECHNICAL NOTE

Detecting Oxidative Activity in Live Cells

Reactive oxygen species (ROS) are formed in cells in response to stimulation and may play a role in signal transduction. The principal probes that have been used to detect ROS in cells have been chemically reduced forms of the xanthenes (fluoresceins and rhodamines) or ethidium, which are colorless and nonfluorescent until oxidized back to the fluorophore by ROS. The precise ROS that mediates the oxidation can be difficult to ascertain. Dichlorodihydrofluorescein diacetate (D-399), dihydroethidium (D-1168) have been the standard dyes for these assays. However, we recommend the chloromethyl derivative of dichlorodihydrofluorescein diacetate (C-6827) because it is better retained in cells before and after stimulation. We also recommend evaluation of our Fc OxyBURST technology for applications that involve following phagocytosis. Reactive OxyBURST reagents permit conjugate of ROS-sensitive probes to other biomolecules to follow their uptake and oxidation by phagocytic cells.




Figure 19.7 Bovine pulmonary artery endothelial cells (BPAEC) were initially stained with the reactive oxygen species (ROS) indicator, 6-carboxy-2',7'-dichlorodihydrofluorescein diacetate, di(acetoxymethyl ester) (C-2938). After a 30-minute incubation, the cells were washed and then incubated simultaneously with FM 5-95 (T-23360) and Hoechst 33342 (H-1399, H-3570, H-21492) in phosphate-buffered saline (PBS) for an additional five minutes before washing and mounting in PBS. The red-fluorescent FM 5-95 appears to stain both the plasma membrane and early endosomes; the green-fluorescent, oxidized carboxydichlorofluorescein localizes to the cytoplasm; and the blue-fluorescent Hoechst 33342 dye stains the nucleus.




Figure 19.8 C-6827 5-(and-6)-chloromethyl-2',7'-dichlorodihydrofluorescein diacetate, acetyl ester.




Figure 19.9 An oxidative burst was detected by flow cytometry of cells labeled with 5-(and-6)-chloromethyl-2',7'-dichlorodihydrofluorescein diacetate, acetyl ester (CM-H₂DCFDA, C-6827). Jurkat cells were incubated with 100 nM CM-H₂DCFDA. The cells were washed and resuspended in either phosphate-buffered saline (PBS, red) or PBS with 0.03% H₂O₂ (blue). The samples were analyzed on a flow cytometer equipped with a 488 nm argon-ion laser and a 525 ± 10 nm bandpass emission filter.

Improved Versions of H₂DCFDA

Intracellular oxidation of H₂DCF tends to be accompanied by leakage of the product, 2',7'-dichlorofluorescein,¹⁹⁴ which may make quantitation or detection of slow oxidation difficult. To enhance retention of the fluorescent product, Molecular Probes offers the carboxylated H₂DCFDA analog (carboxy-H₂DCFDA, C-400), which has two negative charges at physiological pH, and its di(acetoxymethyl ester) (C-2938, Figure 15.79, Figure 19.7), which should more easily pass through membranes during cell loading. Upon oxidation and cleavage of the acetate and ester groups by intracellular esterases, both analogs form carboxydichlorofluorescein (C-368, Section 14.3), with additional negative charges that should impede its leakage out of the cell. Carboxy-H₂DCFDA (C-400) has been used to assess the oxidative process in isolated perfused rat heart tissue¹⁸³ and in transfected *cos-1* cells expressing native or mutagenized prostaglandin endoperoxide H synthase.¹⁹⁵ Its di(acetoxymethyl ester) (C-2938) has been employed to investigate the role of the *bcl-2* proto-oncogene product in preventing apoptosis through its antioxidant properties.¹²⁶ Another probe for oxidative bursts and reactive oxygen species — 5-(and-6)-carboxy-2',7'-difluorodihydrofluorescein diacetate (carboxy-H₂DFFDA, C-13293) — exhibits improved photostability when compared to other fluorescein derivatives in common use.

Molecular Probes has also developed 5-(and-6)-chloromethyl-2',7'-dichlorodihydrofluorescein diacetate, acetyl ester¹⁹⁶ (CM-H₂DCFDA, C-6827; Figure 19.8, Figure 19.9) — a chloromethyl derivative of H₂DCFDA that should exhibit much better retention in live cells. As with our other chloromethyl derivatives (see the description of our Cell-Tracker probes in Section 14.2), we believe that CM-H₂DCFDA passively diffuses into cells, where its acetate groups are cleaved by intracellular esterases and its thiol-reactive chloromethyl group reacts with intracellular glutathione and other thiols. Subsequent oxidation yields a fluorescent adduct that is trapped inside the cell, thus facilitating long-term studies. CM-H₂DCFDA has been used to:

- Measure intracellular reactive oxygen species in cardiac myocytes¹⁹⁷
- Detect insulin-stimulated production of H₂O₂ in insulin-sensitive hepatoma and adipose cells¹⁹⁸
- Detect arsenic-induced production of oxyradicals by confocal laser-scanning microscopy¹⁹⁹
- Explore the parallel induction of reactive oxygen species and Ca²⁺ transients in ouabain-treated myocytes^{197,200}

The diacetate derivatives of the dichlorodihydrofluoresceins are quite stable. When used for intracellular applications, the acetates are cleaved by endogenous esterases, releasing the corresponding dichlorodihydrofluorescein derivative. If, however, these nonfluorescent diacetate derivatives are used for *in vitro* assays, they must first be hydrolyzed with mild base to form the colorless probe.

Dihydrocalcein AM

Calcein AM (C-1430, C-3099, C-3100; Section 15.2, Figure 15.2) is extremely useful as a probe for the study of cell viability, adhesion, multidrug resistance, chemotaxis and other processes. We have combined the superior retention of calcein — the hydrolytic product of calcein AM in viable cells — and the oxidation sensitivity of dihydrofluoresceins to yield the probe dihydrocalcein AM (D-23805, Figure 15.4). Dihydrocalcein AM is freely permeant to cell membranes and is oxidized to green-fluorescent calcein, which has superior retention properties in cells that have intact membranes (Figure 15.3). Dihydrocalcein AM is provided specially packaged as a set of 20 vials, each containing 50 µg of the probe.

OxyBURST Green Reagents

Molecular Probes' Fc OxyBURST Green assay reagent (F-2902) was developed in collaboration with Elizabeth Simons of Boston University to monitor the oxidative burst in phagocytic cells using fluorescence instrumentation. The Fc OxyBURST Green assay reagent comprises bovine serum albumin (BSA) that has been covalently linked to dichlorodihydrofluorescein (H₂DCF) and then complexed with purified rabbit polyclonal anti-BSA antibodies. When these immune complexes bind to Fc receptors, the nonfluorescent H₂DCF molecules are internalized within the phagovacuole and subsequently

oxidized to green-fluorescent dichlorofluorescein (DCF, Figure 16.1; Figure 16.2). Unlike dichlorodihydrofluorescein diacetate (H₂DCFDA), the Fc OxyBURST Green assay reagent does not require intracellular esterases for activation, making this reagent particularly suitable for detecting the oxidative burst in cells with low esterase activity such as monocytes.²⁰¹ The Fc OxyBURST Green assay reagent reportedly produces >8 times more fluorescence than does H₂DCFDA at 60 seconds and >20 times more at 15 minutes following internalization of the immune complex.²⁰²

OxyBURST Green H₂HFF BSA (O-13291) is a sensitive fluorogenic reagent for detecting extracellular release of oxidative products in a spectrofluorometer or a fluorescence microscope (Figure 16.4). This reagent comprises BSA that has been covalently linked to dihydro-2',4,5,6,7,7'-hexafluorofluorescein (H₂HFF), a reduced dye with improved stability. Unlike our Fc OxyBURST Green assay reagent, OxyBURST Green H₂HFF BSA is not complexed with IgG. OxyBURST Green H₂HFF BSA provides up to 1000-fold greater sensitivity than conventional methods based on spectrophotometric detection of superoxide dismutase-inhibitable reduction of cytochrome *c*^{203,204} and allows researchers to take advantage of the sample stirring and temperature control capabilities available in many spectrofluorometers. Because OxyBURST Green H₂HFF BSA is a protein conjugate, it is superior to low molecular weight probes such as dihydrotetramethylrosamine and dihydrorhodamine 123, which are cell permeant and therefore do not exclusively detect extracellular oxidants.

Amine-Reactive OxyBURST Green Reagents

As an alternative to our Fc OxyBURST Green assay reagent and OxyBURST Green H₂HFF BSA, Molecular Probes offers amine-reactive OxyBURST Green H₂DCFDA succinimidyl ester (2',7'-dichlorodihydrofluorescein diacetate, SE; D-2935; Figure 16.6), which can be used to prepare oxidation-sensitive conjugates of a wide variety of biomolecules and particles, including antibodies, antigens, peptides, proteins, dextrans, bacteria, yeast and polystyrene microspheres.^{202,205} Following conjugation to amines, the two acetates of OxyBURST Green H₂DCFDA can be removed by treatment with hydroxylamine at neutral pH to yield the dihydrofluorescein conjugate. The OxyBURST Green H₂DCFDA conjugates are nonfluorescent until they are oxidized to the corresponding fluorescein derivatives. Thus, like our Fc OxyBURST Green assay reagent, they provide a means of detecting the oxidative burst in phagocytic cells. In one application, OxyBURST Green H₂DCFDA succinimidyl ester was conjugated to an antibody that binds specifically to YAC tumor cells. YAC cells opsonized with this customized OxyBURST reagent were then used in a fluorescence microscopy study to show that Fc




Figure 19.10 D-632 dihydrorhodamine 123.

receptor-activated neutrophils appear to deliver reactive oxygen species to the surface of their target cells.²⁰⁵

Dihydrorhodamine 123

Dihydrorhodamine 123 (D-632, D-23806; Figure 19.10, Figure 19.11) is the uncharged and nonfluorescent reduction product of the mitochondrion-selective dye rhodamine 123 (R-302; FluoroPure Grade, R-22420; Section 12.2). This leuco dye passively diffuses across most cell membranes where it is oxidized to cationic rhodamine 123 (Figure 12.23), which localizes in the mitochondria. Like H₂DCF, dihydrorhodamine 123 does not directly detect superoxide,¹⁷² but rather reacts with hydrogen peroxide in the presence of peroxidase,¹⁷² cytochrome *c* or Fe²⁺.¹⁷³ However, dihydrorhodamine 123 also reacts with peroxyxynitrite,^{206,207} the anion formed when nitric oxide reacts with superoxide.^{208,209} Peroxyxynitrite, which may play a role in many pathological conditions,^{62,207} has been shown to react with sulfhydryl groups,²¹⁰ DNA²¹¹ and membrane phospholipids,²¹² as well as with tyrosine²¹³ and other phenolic compounds.²¹⁴

Dihydrorhodamine 123 has been used to investigate reactive oxygen intermediates produced by human and murine phagocytes,²¹⁵ activated rat mast cells²¹⁶ and cultured endothelial cells.¹⁷³ It has also been employed to study the role of the CD14 cell-surface marker in H₂O₂ production by human monocytes.²¹⁷ In addition, dihydrorhodamine 123 has been used with the Fura Red calcium indicator (F-3020, F-3021; Section 20.3) to simultaneously measure oxidative bursts and Ca²⁺ fluxes in monocytes and granulocytes.²¹⁸ Dihydrorhodamine 123 is reportedly a more sensitive probe than H₂DCFDA for detecting granulocyte respiratory bursts.^{219–221}




Figure 19.11 Live bovine pulmonary artery endothelial cells (BPAEC) were first stained with LysoTracker Red DND-99 (L-7528). Then, a solution of dihydrorhodamine 123 (D-632, D-23806) and Hoechst 33258 (H-1398, H-3569, H-21491) was added and allowed to incubate with the cells for an additional 10 minutes before the cells were subsequently washed and visualized. The green-fluorescent oxidation product (rhodamine 123, R-302) localized primarily to the mitochondria. The red-fluorescent LysoTracker Red DND-99 stain accumulated in the lysosomes and the blue-fluorescent Hoechst 33258 dye stained the nuclei. The image was acquired with filters appropriate for DAPI, fluorescein and the Texas Red dye and deconvolved using Huygens software (Scientific Volume Imaging, www.svi.nl).

Dihydrorhodamine 123 is available as a 10 mg vial (D-632) or as a stabilized 5 mM solution in DMSO (D-23806). Because of the susceptibility of dihydrorhodamine 123 to air oxidation, the DMSO solution is recommended when only small quantities are to be used at a time.

A Longer-Wavelength Reduced Rhodamine

Intracellular oxidation of dihydrorhodamine 6G (D-633) yields rhodamine 6G (R-634), which localizes in the mitochondria of living cells (Section 12.2). This cationic oxidation product has longer-wavelength spectra than those of dihydrorhodamine 123, making it especially useful for multicolor applications and in some autofluorescent cells and tissues. Dihydrorhodamine 6G has been used in the study of chloride conductance and mutations of the cystic fibrosis transmembrane conductance regulator.²²²

Reduced MitoTracker Probes

Two of our MitoTracker probes — MitoTracker Orange CM-H₂TMRos (M-7511, Figure 12.10) and MitoTracker Red CM-H₂XRos (M-7513, Figure 12.11) — are chemically reactive reduced rosamines. Unlike MitoTracker Orange CMTMRos and MitoTracker Red CMXRos (M-7510, M-7512; Section 12.2), the reduced versions of these probes do not fluoresce until they enter an actively respiring cell, where they are oxidized by reactive oxygen species to the fluorescent mitochondrion-selective probe and then sequestered in the mitochondria. Both MitoTracker Orange CMTMRos and the reduced MitoTracker Orange CM-H₂TMRos have been used to investigate the metabolic state of *Pneumocystis carinii* mitochondria.²²³




Figure 19.12 Live bovine pulmonary artery endothelial cells (BPAEC) were incubated with the cell-permeant, weakly blue-fluorescent dihydroethidium (D-1168, D-11347, D-23107) and the green-fluorescent mitochondrial stain, MitoTracker Green FM (M-7514). Upon oxidation, red-fluorescent ethidium accumulated in the nucleus.

Dihydroethidium (Hydroethidine)

Although dihydroethidium (Figure 15.20), which is also called hydroethidine, is commonly used to analyze respiratory burst in phagocytes,^{224–226} it has been reported that this probe undergoes significant oxidation in resting leukocytes, possibly through the uncoupling of mitochondrial oxidative phosphorylation.¹⁷¹ Cytosolic dihydroethidium exhibits blue fluorescence; however, once this probe is oxidized to ethidium (E-1305, Section 8.1, Figure 8.12), it intercalates within the cell's DNA, staining its nucleus a bright fluorescent red^{227–229} (Figure 19.12). The mechanism of dihydroethidium's interaction with lysosomes and DNA has been described.²³⁰ Dihydroethidium has been used for many purposes, including:

- Detecting multidrug-resistant cancer cells²³¹ (Section 15.6)
- Following phagocytosis and oxidative bursts by phagocytic blood cells^{224,225,232,233}
- Investigating spermatozoal viability^{234,235}
- Quantitating killer cell–target cell conjugates by flow cytometry methods^{236–238}

Dihydroethidium (hydroethidine) is available in a 25 mg vial (D-1168), as a stabilized 5 mM solution in DMSO (D-23107) or specially packaged in 10 vials of 1 mg each (D-11347); the stabilized DMSO solution or special packaging is recommended when small quantities of the dye will be used over a long period of time.

RedoxSensor Red CC-1 Stain

RedoxSensor Red CC-1 stain (R-14060) is a unique probe whose fluorescence localization appears to be based on a cell's cytosolic redox potential. Scientists at Molecular Probes have found that RedoxSensor Red CC-1 stain passively enters live cells.²³⁹ Once inside, the nonfluorescent probe is either oxidized in the cytosol to a red-fluorescent product (excitation/emission maxima ~540/600 nm), which then accumulates in the mitochondria, or the probe is transported to the lysosomes, where it is oxidized. The differential distribution of the oxidized product between mitochondria and lysosomes appears to depend on the oxidation–reduction (redox) potential of the cytosol. In proliferating cells, mitochondrial staining predominates, whereas contact-inhibited cells primarily exhibit lysosomal staining. The best method we have found to quantitate the distribution of the oxidized product is to use the mitochondrion-selective MitoTracker Green FM stain (M-7514, Section 12.2) in conjunction with the RedoxSensor Red CC-1 stain.

Tetrazolium Salts — Chromogenic Redox Indicators

Tetrazolium salts — especially MTT (M-6494) — are widely used for detecting the redox potential of cells for viability, proliferation and cytotoxicity assays. Upon reduction, these water-soluble colorless compounds form uncharged, brightly colored formazans. Several of the formazans precipitate out of solution and are useful for histochemical localization of the site of reduction or, after solubilization in organic solvent, for quantitation by standard spectrophotometric techniques. The extremely water-soluble formazan product of XTT (X-6493) does not require solubilization prior to quantitation. Many of these salts are re-

duced by specific components of the electron transport chain and may be useful for determining the site of action of specific toxins.²⁴⁰ Selected applications of the tetrazolium salts are listed in Table 19.2. Our Vybrant MTT Cell Proliferation Assay Kit

(V-13154, Section 15.4) provides a means of counting metabolically active cells; the Vybrant MTT assay can detect from 2000 to 250,000 cells, depending on the cell type and conditions. See also Section 15.2 for applications of tetrazolium salts in cells.

References

1. J Photochem Photobiol B 11, 241 (1991);
2. J Am Chem Soc 104, 5541 (1982);
3. J Chem Soc Faraday Trans 1 82, 1627 (1986);
4. Methods Enzymol 186, 635 (1990);
5. Proc Natl Acad Sci USA 76, 6047 (1979);
6. J Am Chem Soc 111, 2909 (1989);
7. Photochem Photobiol 65, 352 (1997);
8. Photochem Photobiol 61, 632 (1995);
9. Photochem Photobiol 61, 529 (1995);
10. Pharmacol Ther 63, 1 (1994);
11. J Photochem Photobiol A 64, 273 (1992);
12. Photochem Photobiol 54, 95 (1991);
13. Photochem Photobiol 53, 169 (1991);
14. Biochem Pharmacol 47, 373 (1994);
15. Biochem Biophys Res Commun 165, 1207 (1989);
16. Biochem Pharmacol 46, 1929 (1993);
17. Antiviral Res 17, 63 (1992);
18. Photochem Photobiol 52, 609 (1990);
19. Photochem Photobiol 43, 677 (1986);
20. J Cell Biol 126, 901 (1994);
21. Adv Photochem 18, 315 (1993);
22. Photochem Photobiol 55, 81 (1992);
23. Proc SPIE-Intl Soc Opt Eng 997, 48 (1988);
24. Photochem Photobiol 66, 374 (1997);
25. Biochim Biophys Acta 1105, 333 (1992);
26. Arch Biochem Biophys 291, 43 (1991);
27. J Photochem Photobiol A 58, 339 (1991);
28. J Photochem Photobiol B 1, 437 (1988);
29. J Immunol Methods 168, 245 (1994);
30. Arch Biochem Biophys 300, 714 (1993);
31. Bone Marrow Transplant 12, 191 (1993);
32. Cancer Res 53, 806 (1993);
33. Free Radic Biol Med 12, 389 (1992);
34. Biochim Biophys Acta 1075, 28 (1991);
35. J Infect Dis 163, 1312 (1991);
36. Photochem Photobiol 53, 1 (1991);
37. Cancer Res 49, 3637 (1989);
38. Biochim Biophys Acta 882, 210 (1986);
39. Methods Enzymol 133, 569 (1986);
40. Biochem Biophys Res Commun 123, 869 (1984);
41. Arch Biochem Biophys 314, 284 (1994);
42. J Reprod Fertil 97, 441 (1993);
43. Meth Enzymol 12, 5 (1967);
44. Anal Biochem 290, 138 (2001);
45. J Biol Chem 273, 32030 (1998);
46. Curr Biol 6, 1497 (1996);
47. Methods Cell Biol 44, 715 (1994);
48. Proc Natl Acad Sci U S A 91, 2659 (1994);
49. Nucleic Acids Res 22, 4789 (1994);
50. Bioconjug Chem 4, 69 (1993);
51. Methods Enzymol 208, 414 (1991);
52. Annu Rev Biochem 59, 207 (1990);
53. J Am Chem Soc 111, 4941 (1989);
54. Proc Natl Acad Sci U S A 86, 9702 (1989);
55. Biochemistry 29, 8447 (1990);
56. J Am Chem Soc 109, 1990 (1987);
57. Anal Biochem 212, 85 (1993);
58. Anal Chem 68, 867 (1996);
59. FASEB J 9, 1085 (1995);
60. Anal Chem 69, 4295 (1997);
61. J Biol Chem 276, 35253 (2001);
62. Proc Natl Acad Sci U S A 87, 1620 (1990);
63. Eur J Biochem 221, 695 (1994);
64. Anal Biochem 206, 273 (1992);
65. Arch Biochem Biophys 373, 447 (2000);
66. Circ Res 84, 1203 (1999);
67. Biochem Biophys Res Commun 248, 382 (1998);
68. J Biol Chem 273, 33972 (1998);
69. Inflammation 20, 151 (1996);
70. J Biolumin Chemilumin 12, 277 (1997);
71. J Immunol Methods 97, 209 (1987);
72. Blood 83, 3324 (1994);
73. Free Radic Biol Med 15, 447 (1993);
74. Cytometry 18, 147 (1994);
75. Free Radic Biol Med 17, 117 (1994);
76. J Cell Physiol 151, 466 (1992);
77. Biochem J 289, 587 (1993);
78. Free Radic Biol Med 15, 143 (1993);
79. Am J Physiol 266, G799 (1994);
80. Life Sci 55, 1427 (1994);
81. Chem Soc Rev 6, 195 (1977);
82. Can J Biochem 51, 158 (1973);
83. Biochim Biophys Acta 385, 232 (1975);
84. Biochem Biophys Res Commun 163, 836 (1989);
85. Biochim Biophys Acta 981, 235 (1989);
86. FEBS Lett 169, 169 (1984);
87. J Cell Biol 273, 26900 (1998);
88. Neurochem Res 22, 333 (1997);
89. Am J Physiol 271, F209 (1996);
90. Biochim Biophys Acta 1487, 61 (2000);
91. Circ Res 84, 229 (1999);
92. J Biol Chem 276, 15575 (2001);
93. Plant Physiol 125, 1591 (2001);
94. Atherosclerosis 152, 307 (2000);
95. J Biol Chem 269, 798 (1994);
96. J Biol Chem 274, 2234 (1999);
97. Free Radic Biol Med 11, 81 (1991);
98. J Exp Med 139, 208 (1974);
99. Free Radic Biol Med 22, 93 (1997);
100. Proc Soc Exp Biol Med 206, 53 (1994);
101. Biochem Int 10, 205 (1985);
102. J Invest Dermatol 110, 966 (1998);
103. J Immunol Methods 202, 133 (1997);
104. FASEB J 13, 69 (1999);
105. Proc Natl Acad Sci U S A 97, 10930 (2000);
106. Science 293, 1806 (2001);
107. Anal Biochem 253, 162 (1997);
108. J Neurochem 79, 266 (2001);
109. Am J Physiol Lung Cell Mol Physiol 281, L993 (2001);
110. Mol Hum Reprod 7, 237 (2001);
111. J Invest Dermatol 112, 751 (1999);
112. Clin Chim Acta 63, 37 (1975);
113. Free Radic Biol Med 8, 121 (1990);
114. J Biol Chem 273, 7828 (1998);
115. Biochemistry 37, 13781 (1998);
116. Biochim Biophys Acta 1330, 127 (1997);
117. J Biochem Biophys Methods 35, 23 (1997);
118. J Biol Chem 272, 12328 (1997);
119. Biochemistry 34, 12755 (1995);
120. Proc Natl Acad Sci U S A 91, 1183 (1994);
121. Arch Biochem Biophys 297, 147 (1992);
122. Cytometry 13, 686 (1992);
123. Anal Biochem 196, 443 (1991);
124. Chem Phys Lipids 53, 309 (1990);
125. Biochim Biophys Acta 944, 29 (1988);
126. Cell 75, 241 (1993);
127. Biochemistry 16, 819 (1977);
128. Proc Natl Acad Sci USA 77, 26 (1980);
129. Biosci Biotechnol Biochem 56, 605 (1992);
130. Anal Lett 21, 965 (1988);
131. J Chromatogr A 881, 159 (2000);
132. Lipids 31, 1091 (1996);
133. J Chromatogr 617, 205 (1993);
134. J Chromatogr 622, 153 (1993);
135. FEBS Lett 474, 137 (2000);
136. Free Radic Res 29, 43 (1998);
137. J Agric Food Chem 48, 1150 (2000);
138. Anal Biochem 265, 290 (1998);
139. FEBS Lett 453, 278 (1999);
140. Biochim Biophys Acta 1298, 119 (1996);
141. Free Radic Biol Med 21, 241 (1996);
142. Biochem Pharmacol 49, 1649 (1995);
143. Clin Chem 41, 1738 (1995);
144. Biochemistry 31, 8090 (1992);
145. Anal Biochem 177, 300 (1989);
146. Life Sci 60, PL 301 (1997);
147. Free Radic Biol Med 18, 29 (1995);
148. J Biol Chem 273, 33279 (1998);
149. FASEB J 2, 2487 (1988);
150. Biochem Mol Biol Int 33, 1179 (1994);
151. Free Radic Biol Med 6, 623 (1989);
152. J Immunol 129, 1589 (1982);
153. J Appl Physiol 76, 539 (1994);
154. J Leukoc Biol 54, 300 (1993);
155. J Biochem (Tokyo) 106, 355 (1989);
156. Biochem Biophys Res Commun 155, 106 (1988);
157. Am J Physiol 262, G1015 (1992);
158. Arch Biochem Biophys 310, 352 (1994);
159. Anal Chem 65, 1794 (1993);
160. Biochem J 290, 51 (1993);
161. Free Radic Biol Med 18, 1 (1995);
162. J Lipid Res 33, 1051 (1992);
163. Methods Enzymol 233, 324 (1994);
164. Lipids 33, 1235 (1998);
165. Biochim Biophys Acta 1097, 145 (1991);
166. Biomed Biochim Acta 49, 991 (1990);
167. Proc Natl Acad Sci U S A 89, 11426 (1992);
168. Arch Toxicol 68, 582 (1994);
169. Brain Res 635, 113 (1994);
170. Chem Res Toxicol 5, 227 (1992);
171. J Leukoc Biol 47, 440 (1990);
172. Eur J Biochem 217, 973 (1993);
173. Arch Biochem Biophys 302, 348 (1993);
174. Exp Cell Res 209, 375 (1993);
175. J Immunol Methods 159, 173 (1993);
176. J Immunol Methods 159, 131 (1993);
177. Cytometry 13, 615 (1992);
178. Cytometry 13, 525 (1992);
179. Methods Enzymol 233, 128 (1994);
180. Lab Invest 64, 167 (1991);
181. Am J Physiol 264, H881 (1993);
182. Lab Invest 70, 579 (1994);
183. Free Radic Res Commun 16, 217 (1992);
184. Cytometry 14, 287 (1993);
185. Cell 77, 817 (1994);
186. Brain Res 610, 172 (1993);
187. J Neurosci 14, 2260 (1994);
188. Eur J Cell Biol 65, 200 (1994);
189. J Biol Chem 269, 14546 (1994);
190. Toxicol Appl Pharmacol 120, 29 (1993);
191. J Cell Biol 126, 1079 (1994);
192. J Leukoc Biol 51, 496 (1992);
193. J Biol Chem 274, 28161 (1999);
194. Free Radic Biol Med 16, 509 (1994);
195. Biochemistry 34, 7194 (1995);
196. We have discovered by mass spectroscopy that the current lot and all previous lots of this product unexpectedly contain an extra acetate moiety that is covalently bound to the carboxylic acid. Consequently the product's name, molecular weight and chemical structure have been changed to reflect presence of this residue. This change does not affect the use of this material for its intended purposes and the product is in every way identical to material previously provided by Molecular Probes, which is the exclusive supplier of this material.
197. J Biol Chem 274, 19323 (1999);
198. J Biol Chem 276, 21938 (2001);
199. Proc Natl Acad Sci U S A 98, 1643 (2001);
200. J Biol Chem 275, 27838 (2000);
201. J Leukoc Biol 43, 304 (1988);
202. J Immunol Methods 130, 223 (1990);
203. J Biol Chem 255, 1874 (1980);
204. J Clin Invest 61, 1081 (1978);
205. J Cell Physiol 156, 428 (1993);
206. Biochemistry 34, 3544 (1995);
207. Free Radic Biol Med 16, 149 (1994);
208. Free Radic Res Commun 18, 195 (1993);
209. Free Radic

References — continued

- Res Commun 10, 221 (1990); **210.** J Biol Chem 266, 4244 (1991); **211.** J Am Chem Soc 114, 5430 (1992); **212.** Arch Biochem Biophys 288, 481 (1991); **213.** Arch Biochem Biophys 298, 431 (1992); **214.** Arch Biochem Biophys 298, 438 (1992); **215.** J Immunol Methods 131, 269 (1990); **216.** APMIS 102, 474 (1994); **217.** FEBS Lett 273, 55 (1990); **218.** Cytometry 13, 693 (1992); **219.** J Immunol Methods 178, 89 (1995); **220.** J Immunol Methods 162, 261 (1993); **221.** Naturwissenschaften 75, 354 (1988); **222.** Proc Natl Acad Sci U S A 93, 1167 (1996); **223.** J Euk Microbiol 41, 79 (1994); **224.** FEMS Microbiol Lett 122, 187 (1994); **225.** J Immunol Methods 170, 117 (1994); **226.** Cytometry 12, 687 (1991); **227.** J Leukoc Biol 55, 253 (1994); **228.** J Histochem Cytochem 34, 1109 (1986); **229.** Biotechniques 3, 270 (1985); **230.** Histochemistry 94, 205 (1990); **231.** Exp Cell Res 190, 69 (1990); **232.** Cytometry 17, 294 (1994); **233.** J Immunol 151, 1463 (1993); **234.** J Histochem Cytochem 39, 485 (1991); **235.** Gamete Res 22, 355 (1989); **236.** Cytometry 12, 666 (1991); **237.** J Immunol Methods 139, 281 (1991); **238.** J Immunol Methods 130, 251 (1990); **239.** Free Radic Biol Med 28, 1266 (2000); **240.** Toxicity Testing Using Microorganisms, Vol. I, Bitton G, Dutka BJ, Eds. pp. 27–55 (1986).

Data Table — 19.2 Generating and Detecting Reactive Oxygen Species

| Cat # | MW | Storage | Soluble | Abs | EC | Em | Solvent | Notes |
|---------|-----------|------------|-------------------------|------|-----------|-----------|-------------------|------------|
| A-7923 | 376.48 | F,D,L | DMSO | 358 | 11,000 | 424 | MeOH | 1 |
| A-12222 | 257.25 | FF,D,A | DMSO | 280 | 6,000 | none | pH 8 | 2 |
| B-3932 | 448.32 | F,L | DMSO, CHCl ₃ | 665 | 161,000 | 676 | MeOH | |
| C-400 | 531.30 | F,D | DMSO, EtOH | 290 | 5,600 | none | MeCN | 3 |
| C-2938 | 675.43 | F,D,AA | DMSO | 291 | 5,700 | none | MeOH | 3 |
| C-2944 | 423.47 | FF,D,LL,AA | MeOH | 429 | 7,500 | see Notes | pH 7 | 4, 5, 6 |
| C-6827 | 577.80 | F,D,AA | DMSO | 287 | 9,100 | none | MeOH | 3 |
| C-7924 | 487.62 | F,D,L | DMSO, H ₂ O | 385 | 5,800 | 485 | pH 7 | 1 |
| C-13293 | 498.39 | F,D | DMSO, EtOH | 290 | 5,500 | none | MeCN | 7 |
| D-399 | 487.29 | F,D | DMSO, EtOH | 258 | 11,000 | none | MeOH | 3 |
| D-632 | 346.38 | F,D,L,AA | DMF, DMSO | 289 | 7,100 | none | MeOH | 8, 9 |
| D-633 | 444.57 | F,D,L,AA | DMF, DMSO | 296 | 11,000 | none | MeOH | 8, 9 |
| D-1168 | 315.42 | FF,L,AA | DMF, DMSO | 355 | 14,000 | see Notes | MeCN | 8, 10 |
| D-2935 | 584.37 | F,D,AA | DMF | 258 | 11,000 | none | MeOH | 3 |
| D-3861 | 504.43 | F,L | DMSO | 582 | 140,000 | 591 | MeOH | 11 |
| D-7894 | 386.43 | F,D,LL | MeCN | 358 | 29,000 | none | MeOH | 12 |
| D-11347 | 315.42 | FF,L,AA | DMF, DMSO | 355 | 14,000 | see Notes | MeCN | 8, 10 |
| D-23107 | 315.42 | FF,D,L,AA | DMSO | 355 | 14,000 | see Notes | MeCN | 10, 13 |
| D-23805 | 1068.95 | F,D | DMSO | 285 | 5,800 | none | MeCN | 14 |
| D-23806 | 346.38 | F,D,L,AA | DMSO | 289 | 7,100 | none | MeOH | 9, 13 |
| F-2902 | see Notes | RR,L,AA | H ₂ O | <300 | | none | | 15, 16, 17 |
| H-7476 | 504.45 | F,D,L | DMSO, DMF | 591 | 37,000 | 594 | EtOH | |
| H-7515 | 546.53 | F,D,L | DMSO | 463 | 24,000 | 600 | MeOH | 18 |
| H-7516 | 528.51 | F,D,L | DMSO | 459 | 23,000 | 616 | MeOH | 18 |
| L-6868 | 510.50 | L | H ₂ O | 455 | 7,400 | 505 | H ₂ O | 19, 20 |
| L-8455 | 177.16 | D,L | DMF | 355 | 7,500 | 411 | MeOH | 20 |
| M-689 | 485.98 | F,DD,L | DMF, DMSO | 629 | 75,000 | none | MeCN | 21 |
| M-6494 | 414.32 | D,L | H ₂ O, DMSO | 375 | 8,300 | none | MeOH | 22, 23 |
| M-7511 | 392.93 | F,D,L,AA | DMSO | 235 | 57,000 | none | MeOH | 8, 9 |
| M-7513 | 497.08 | F,D,L,AA | DMSO | 245 | 45,000 | none | MeOH | 8, 9 |
| M-7913 | 258.32 | F,L | DMF, DMSO | 352 | 30,000 | 401 | MeOH | 24 |
| M-23800 | 291.74 | FF,D,LL,AA | DMSO | 430 | 8,400 | 546 | MeOH | 25 |
| M-24571 | 569.67 | D,L | DMSO, EtOH | 555 | 143,000 | 578 | MeOH | |
| N-6495 | 817.65 | D,L | H ₂ O, DMSO | 256 | 64,000 | none | MeOH | 22 |
| P-800 | ~240,000 | RR,L | see Notes | 546 | 2,410,000 | 575 | pH 7 | 26 |
| P-801 | ~240,000 | RR,L | see Notes | 565 | 1,960,000 | 578 | pH 7 | 26 |
| P-1901 | 276.42 | FF,LL,AA | EtOH | 303 | 76,000 | 416 | MeOH | 27 |
| P-6879 | 363.16 | F,D,L | DMSO | 270 | 28,000 | none | CHCl ₃ | 28 |
| R-14000 | 1057.75 | F,D | DMSO | 313 | 9,700 | none | MeOH | 29 |
| R-14060 | 434.41 | F,D,L,AA | DMSO | 239 | 52,000 | none | MeOH | 8, 30 |
| T-7539 | 740.87 | F,D,L | MeOH, DMF | 524 | 91,000 | 550 | MeOH | |
| X-6493 | 674.53 | D | H ₂ O, DMSO | 286 | 15,000 | none | MeOH | 31 |

For definitions of the contents of this data table, see "How to Use This Book" on page viii.

Notes

1. Fluorescence of A-7923 and C-7924 is weak. Reaction of the nitroxide moiety with superoxide or hydroxyl radicals results in increased fluorescence without a spectral shift (Anal Biochem 212, 85 (1993)).
2. Peroxidase-catalyzed reaction of A-12222 with H₂O₂ produces fluorescent resorufin R-363 (Section 10.1).
3. Dihydrofluorescein diacetates are colorless and nonfluorescent until both the acetates are hydrolyzed and the products are subsequently oxidized to fluorescein derivatives. The materials contain less than 0.1% of oxidized derivative when initially prepared. The end products from C-400, C-2938, C-6827, D-399 and D-2935 are 2',7'-dichlorofluorescein derivatives with spectra similar to C-368 (Section 21.3).
4. C-2944 emits chemiluminescence (Em = 466 nm) on oxidation by superoxide (Anal Biochem 206, 273 (1992)).
5. Do NOT dissolve in DMSO.
6. Aqueous solutions of coelenterazine (>1 mM) can be prepared in pH 7 buffer containing 50 mM 2-hydroxypropyl-β-cyclodextrin (Biosci Biotechnol Biochem 61, 1219 (1997)).
7. Difluorodihydrofluorescein diacetates are colorless and nonfluorescent. Acetate hydrolysis and subsequent oxidation generate a fluorescent difluorofluorescein with spectra similar to O-6146 (Section 1.5).

8. This compound is susceptible to oxidation, especially in solution. Store solutions under argon or nitrogen. Oxidation appears to be catalyzed by illumination.
9. These compounds are essentially colorless and nonfluorescent until oxidized. Oxidation products (in parentheses) are as follows: D-632 and D-23806 (R-302; Section 12.2); D-633 (R-634; Section 12.2); M-7511 (M-7510; Section 12.2); M-7513 (M-7512; Section 12.2).
10. Dihydroethidium has blue fluorescence (Em ~420 nm) until oxidized to ethidium E-1305 (Section 8.1). The reduced dye does not bind to nucleic acids (FEBS Lett 26, 169 (1972)).
11. Oxidation of the polyunsaturated butadienyl portion of the BODIPY 581/591 dye results in a shift of the fluorescence emission peak from ~590 nm to ~510 nm (Methods Enzymol 319, 603 (2000); FEBS Lett 453, 278 (1999)).
12. Oxidation product is strongly fluorescent. Em = 379 nm. Oxidation occurs rapidly in solution when illuminated.
13. This product is supplied as a ready-made solution in DMSO with sodium borohydride added to inhibit oxidation.
14. D-23805 is colorless and nonfluorescent until the AM ester groups are hydrolyzed and the resulting leuco dye is subsequently oxidized. The final product is calcein (C-481, Section 14.3).
15. This product is supplied as a ready-made solution in the solvent indicated under **Soluble**.
16. F-2902 is essentially colorless and nonfluorescent until oxidized. A small amount (~5%) of oxidized material is normal and acceptable for the product as supplied. The oxidation product is fluorescent (Abs = 495 nm, Em = 524 nm) (J Immunol Methods 130, 223 (1990)).
17. This product consists of a dye-bovine serum albumin conjugate (MW ~66,000) complexed with IgG in a ratio of approximately 1:4 mol:mol (BSA:IgG).
18. H-7515 and H-7516 have weaker absorption peaks at longer wavelengths: Abs = 581 nm (EC = 12,000 cm⁻¹M⁻¹) for H-7515, Abs = 580 nm (EC = 9000 cm⁻¹M⁻¹) for H-7516. Their quantum yields for singlet oxygen production (0.83 for H-7515, 0.76 for H-7516) are essentially independent of the irradiation wavelength between 430 and 580 nm (J Photochem Photobiol A 64, 273 (1992)).
19. L-6868 has much stronger absorption at shorter wavelengths (Abs = 368 nm (EC = 36,000 cm⁻¹M⁻¹)).
20. This compound emits chemiluminescence upon oxidation in basic aqueous solutions. Emission peaks are at 425 nm (L-8455) and 470 nm (L-6868).
21. Isothiocyanates are unstable in water and should not be stored in aqueous solution.
22. Enzymatic reduction products are water-insoluble formazans with Abs = 505 nm (M-6494) and 605 nm (N-6495) after solubilization in DMSO or DMF. See literature sources for further information (Histochemistry 76, 381 (1982); Prog Histochem Cytochem 9, 1 (1976)).
23. M-6494 also has Abs = 242 nm (EC = 21,000 cm⁻¹M⁻¹) in MeOH.
24. Generates chemiluminescence (Em = 465 nm in 0.1 M SDS) upon reaction with ¹O₂ (J Am Chem Soc 108, 4498 (1986)).
25. Generates chemiluminescence (Em = 455 nm) upon reaction with superoxide.
26. Phycobiliproteins are packaged as suspensions in 60% ammonium sulfate, pH 7.0. Store refrigerated at 4°C but DO NOT FREEZE.
27. *cis*-Parinaric acid is readily oxidized to nonfluorescent products. Use under N₂ or Ar except when oxidation is intended. Stock solutions should be prepared in deoxygenated solvents. *cis*-Parinaric acid is appreciably fluorescent in lipid environments and organic solvents but is nonfluorescent in water.
28. Iodoacetamides in solution undergo rapid photodecomposition to unreactive products. Minimize exposure to light prior to reaction.
29. Acetate hydrolysis of R-14000 yields rose bengal (Abs = 556 nm (EC = 104,000 cm⁻¹M⁻¹) Em = 572 nm in MeOH) (Photochem Photobiol 66, 374 (1997)).
30. R-14060 is colorless and nonfluorescent until oxidized. The spectral characteristics of the oxidation product (2,3,4,5,6-pentafluorotetramethylrosamine) are similar to those of T-639 (Section 12.2).
31. Enzymatic reduction product is a water-soluble formazan, Abs = 475 nm.


Product List — 19.2 Generating and Detecting Reactive Oxygen Species

| Cat # | Product Name | Unit Size |
|---------|---|------------|
| A-7923 | 4-((9-acridinecarbonyl)amino)-2,2,6,6-tetramethylpiperidin-1-oxyl, free radical (TEMPO-9-AC) | 5 mg |
| A-22177 | Amplex [®] Red reagent *packaged for high-throughput screening* | 10 x 10 mg |
| A-22188 | Amplex [®] Red Hydrogen Peroxide/Peroxidase Assay Kit *500 assays* | 1 kit |
| A-12222 | Amplex [®] Red reagent (10-acetyl-3,7-dihydroxyphenoxazine) | 5 mg |
| A-22182 | Amplex [®] Red Xanthine/Xanthine Oxidase Assay Kit *400 assays* | 1 kit |
| B-3932 | (<i>E,E</i>)-3,5-bis-(4-phenyl-1,3-butadienyl)-4,4-difluoro-4-bora-3a,4a-diaza-s-indacene (BODIPY [®] 665/676) | 5 mg |
| C-400 | 5-(and-6)-carboxy-2',7'-dichlorodihydrofluorescein diacetate (carboxy-H ₂ DCFDA) *mixed isomers* | 25 mg |
| C-2938 | 6-carboxy-2',7'-dichlorodihydrofluorescein diacetate, di(acetoxymethyl ester) | 5 mg |
| C-13293 | 5-(and-6)-carboxy-2',7'-difluorodihydrofluorescein diacetate (carboxy-H ₂ DFFDA) *mixed isomers* | 5 mg |
| C-7924 | 5-(2-carboxyphenyl)-5-hydroxy-1-((2,2,5,5-tetramethyl-1-oxypyrrolidin-3-yl)methyl)-3-phenyl-2-pyrrolin-4-one, potassium salt (proxyl fluorescamine) | 5 mg |
| C-6827 | 5-(and-6)-chloromethyl-2',7'-dichlorodihydrofluorescein diacetate, acetyl ester (CM-H ₂ DCFDA) *mixed isomers* *special packaging* | 20 x 50 µg |
| C-2944 | coelenterazine | 250 µg |
| D-399 | 2',7'-dichlorodihydrofluorescein diacetate (2',7'-dichlorofluorescein diacetate; H ₂ DCFDA) | 100 mg |
| D-2935 | 2',7'-dichlorodihydrofluorescein diacetate, succinimidyl ester (OxyBURST [®] Green H ₂ DCFDA, SE) | 5 mg |
| D-3861 | 4,4-difluoro-5-(4-phenyl-1,3-butadienyl)-4-bora-3a,4a-diaza-s-indacene-3-undecanoic acid (BODIPY [®] 581/591 C ₁₁) | 1 mg |
| D-23805 | dihydrocalcein, AM *special packaging* | 20 x 50 µg |
| D-1168 | dihydroethidium (hydroethidine) | 25 mg |
| D-11347 | dihydroethidium (hydroethidine) *special packaging* | 10 x 1 mg |
| D-23107 | dihydroethidium (hydroethidine) *5 mM stabilized solution in DMSO* | 1 mL |
| D-632 | dihydrorhodamine 123 | 10 mg |
| D-23806 | dihydrorhodamine 123 *5 mM stabilized solution in DMSO* | 1 mL |
| D-633 | dihydrorhodamine 6G | 25 mg |
| D-7894 | diphenyl-1-pyrenylphosphine (DPPP) | 5 mg |
| F-2902 | Fc OxyBURST [®] Green assay reagent *25 assays* *3 mg/mL* | 500 µL |
| H-7476 | hypericin | 1 mg |
| H-7515 | hypocrellin A | 1 mg |
| H-7516 | hypocrellin B | 1 mg |
| L-6868 | lucigenin (bis- <i>N</i> -methylacridinium nitrate) *high purity* | 10 mg |
| L-8455 | luminol (3-aminophthalhydrazide) | 25 g |
| M-689 | malachite green isothiocyanate | 10 mg |

| Cat # | Product Name | Unit Size |
|---------|---|------------|
| M-24571 | merocyanine 540 | 25 mg |
| M-7913 | <i>trans</i> -1-(2'-methoxyvinyl)pyrene | 1 mg |
| M-23800 | 2-methyl-6-(4-methoxyphenyl)-3,7-dihydroimidazo[1,2-a]pyrazin-3-one, hydrochloride (MCLA) | 5 mg |
| M-7511 | MitoTracker [®] Orange CM-H ₂ TMRos *special packaging* | 20 x 50 µg |
| M-7513 | MitoTracker [®] Red CM-H ₂ XRos *special packaging* | 20 x 50 µg |
| M-6494 | MTT (3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide) | 1 g |
| N-6495 | nitro blue tetrazolium chloride (NBT) | 1 g |
| O-13291 | OxyBURST [®] Green H ₂ HFF BSA *special packaging* | 5 x 1 mg |
| P-1901 | <i>cis</i> -parinaric acid *special packaging* | 10 x 10 mg |
| P-6879 | <i>N</i> -(1,10-phenanthrolin-5-yl)iodoacetamide | 5 mg |
| P-800 | B-phycoerythrin *4 mg/mL* | 0.5 mL |
| P-801 | R-phycoerythrin *4 mg/mL* | 0.5 mL |
| R-14060 | RedoxSensor [™] Red CC-1 *special packaging* | 10 x 50 µg |
| R-14000 | rose bengal diacetate | 5 mg |
| T-7539 | 2',4',5',7'-tetrabromorhodamine 123 bromide | 5 mg |
| X-6493 | XTT (2,3-bis-(2-methoxy-4-nitro-5-sulfophenyl)-2 <i>H</i> -tetrazolium-5-carboxanilide) | 100 mg |

19.3 Probes for Nitric Oxide Research

Nitric oxide (NO), the molecule that makes the firefly glow,^{1,2} plays a critical role as a molecular mediator of a variety of physiological processes, including blood-pressure regulation and neurotransmission.³⁻⁸ In endothelial cells, as well as neurons and astrocytes, NO is synthesized from L-arginine in a reaction catalyzed by nitric-oxide synthase (NOS) (Figure 19.13). NO that diffuses into smooth muscle cells binds to the heme group of guanylate cyclase. Because free NO is a transient species with a half-life of about five seconds, many investigations of this gaseous molecule have relied largely on studies of NOS. Preparing NO solutions and detecting NO in experimental systems require special precautions to achieve reproducibility.⁹ NO can also complex with superoxide to form the strong oxidant, peroxy nitrite anion¹⁰ (ONOO⁻, Table 19.1), which may be a major cytotoxic agent produced during inflammation, sepsis and ischemia/reperfusion.¹¹ Activated macrophage and neutrophils produce nitric oxide and superoxide, and thus peroxy nitrite anion, at similar rates.¹² NO generators are also reported to produce an accumulation of chelatable Zn²⁺ in hippocampal neuronal perikarya, as determined with some of our Zn²⁺ indicators¹³ (Section 20.7, Table 20.6).



Spontaneous Nitric Oxide Donors and Antagonist

DEANO and Spermine NONOate

The DEANO (diethylamine nitric oxide, D-7915) and spermine NONOate (S-7916) solids provide a means of preparing aqueous NO solutions.¹⁴ When these solids are dissolved in buffer, cell culture medium or blood, they dissociate to form two molecules of NO and one molecule of the corresponding amine¹⁵ (Figure 19.14). The delivery of NO can be easily controlled by preparing moderately basic solutions of these NONOates and then lowering the pH to initiate NO generation.

DEANO has a half-life of about two minutes in pH 7.4 phosphate buffer at 37°C, releasing two molecules of NO and one molecule of diethylamine.¹⁴⁻¹⁷ DEANO therefore allows more sustained delivery of NO than is possible using NO-saturated aqueous solutions (Figure 19.15). The data in Figure 19.15 are from a study of cytochrome P-450 inhibition by NO in which DEANO was used in conjunction with fluorogenic microsomal dealkylase substrates¹⁸ (R-441, R-352; Section 10.6). DEANO has also been used as a NO donor to stimulate ADP-ribosylation of various proteins.¹⁹

With a half-life of 39 minutes at 37°C in pH 7.4 buffer, spermine NONOate releases NO much more slowly than does DEANO, making it suitable for whole animal infusions and experiments with long incubations.¹⁵ Also, spermine, the by-product of spermine NONOate dissociation, is less likely to be deleterious in living systems and may also demonstrate biological activity of its own.²⁰ Both DEANO and spermine NONOate are reported to inhibit *in vitro* proliferation of A375 melanoma cells.²¹

SNAP and SIN-1

Molecular Probes also offers the NO donors SNAP (*S*-nitroso-*N*-acetylpenicillamine, N-7892, N-7927) and SIN-1 (3-morpholininosydnimine, hydrochloride; M-7891, M-7926), which spontaneously release NO (and superoxide in the case of SIN-1) under